



ARROW PHARMATECH
EXPERTISE . SERVICE . INTEGRITY

ARROW PHARMATECH



BULK DENSITY TESTER VOLUMETER (BT-200i) Premium (iSmart)

SPECIFICATIONS

MODEL	BT-200i (-iSmart)
TOP FUNNEL WITH MESH	18 mesh
DISPLAY	LCD 20 X 4 character
BAFFLE / MATERIAL	4 NOS / GLASS
BAFFLE BOX	45mm X 45 mm
BOTTOM FUNNEL	SS 316
INSTRUMENT STAND	LEVEL & VIBRATION FREE
CYLINDRICAL CUP	SS316 Cylindrical receiving cup with an volume of 25.00 ± 0.05 ml and inside diameter of 30.0±2.0mm
LOAD CELL	5 KG ACCURACY 0.01 mg
Power Supply	230 V AC, 50 Hz, 880 VA

ARROW PHARMATECH

C-004, Vijaya CHS., D.V. Marg, Vitawa, Thane (W) - 400 605. India

Email : arrowpharmatech@gmail.com / Contact No. 8169231891 / 7400063721

URL: www.arrowpharmatech.com E-mail: arrowpharmatech@gmail.com, kumar.m@arrowpharmatech.com