



ARROW PHARMATECH
EXPERTISE . SERVICE . INTEGRITY

ARROW PHARMATECH



40 Mesh Basket



10 Mesh Basket



20 Mesh Basket



100 Mesh Basket



150 Mesh Basket



Suppository Basket



Stationary Basket



Alfuzosin Tablet Holder



Dissolution Paddles



'Sotax Style' Spiral Sinkers



Basket Sinkers



Basket Film Sinkers

ARROW PHARMATECH

C-004, Vijaya CHS., D.V. Marg, Vitawa, Thane (W) - 400 605. India

Email : arrowpharmatech@gmail.com / Contact No. 8169231891 / 7400063721

URL: www.arrowpharmatech.com E-mail: arrowpharmatech@gmail.com, kumar.m@arrowpharmatech.com



ARROW PHARMATECH

ARROW PHARMATECH
EXPERTISE . SERVICE . INTEGRITY



Japanese Sinker
(Figure 2a Ref USP Above)



Sinker or Coil Sinker



Capsule Sinker
STD



3 Prong Sinker



Sinkers



Metformin Sinkers



Paddle Over Disc



3 Prong Sinker

Capsule Size Guide

[← Back to Dissolution Accessories](#)

Capsule Size Guide

| Capsule Size | Length | | Diameter | |
|--------------|----------|------|----------|------|
| | (inches) | (mm) | (inches) | (mm) |
| 000 | 1.029 | 26.1 | 0.390 | 9.91 |
| 00E | 1.012 | 25.7 | 0.336 | 8.53 |
| 00 | 0.917 | 23.3 | 0.336 | 8.53 |
| 0E | 0.909 | 23.1 | 0.301 | 7.65 |
| 0 | 0.854 | 21.7 | 0.301 | 7.65 |
| 1 | 0.765 | 19.4 | 0.272 | 6.91 |
| 2 | 0.709 | 18.0 | 0.250 | 6.35 |
| 3 | 0.626 | 15.9 | 0.229 | 5.83 |
| 4 | 0.563 | 14.3 | 0.209 | 5.32 |
| 5 | 0.437 | 11.1 | 0.193 | 4.91 |

ARROW PHARMATECH

C-004, Vijaya CHS., D.V. Marg, Vitawa, Thane (W) - 400 605. India

Email : arrowpharmatech@gmail.com / Contact No. 8169231891 / 7400063721

URL: www.arrowpharmatech.com E-mail: arrowpharmatech@gmail.com, kumar.m@arrowpharmatech.com



ARROW PHARMATECH

ARROW PHARMATECH
EXPERTISE . SERVICE . INTEGRITY



Pharmed Tube for Ismatech Pump



Tygon Tube for Ismatech Pump



Vials



Vials



SS Filter




10 Micron Filter With Holder



Bowl Stand 1



Bowl Stand 2



250 ml Dissolution Vessel



2000 ml Amber Vessel



1 Liter Clear



Peak 1



Amber 1



Amber 2



Amber 3



2LITRE



Diffusion Cell



DT Basket Assembly



HTC Vibration Meter

ARROW PHARMATECH

C-004, Vijaya CHS., D.V. Marg, Vitawa, Thane (W) - 400 605. India

Email : arrowpharmatech@gmail.com / Contact No. 8169231891 / 7400063721

URL: www.arrowpharmatech.com E-mail: arrowpharmatech@gmail.com, kumar.m@arrowpharmatech.com