



ARROW PHARMATECH
EXPERTISE . SERVICE . INTEGRITY

ARROW PHARMATECH

Electromagnetic Sieve Shakers

SOP STORAGE UPTO 16 PROGRAMS

SIMULTANEOUS MODE PROGRAM OPTION

DIGITAL AMPLITUDE



10mm
20µm

NEW



Premium Model (iPremium-QS-200P)

Dry Sieving	Standard
Wet Sieving	Optional Kit
Mode of Operation	3 modes (Continuous, Intermittent & Advance)
Intermittent Operation	At intervals upto 0.5 seconds
Advance Mode	Programmable time for Continuous & Intermittent (Together) in this mode
Capacity	Upto 8 Sieves of 50mm height Upto 16sieves of 25 mm height
Shake Time	Programmable from 1 mm to 99mm
Digital Amplitude Display	Programmable from 0.1mm to 3.0mm 20 x 4 character alphanumeric LCD
Noise Level	Less than 61 db without sieves at maximum amplitude. Less than 71 db with sieves and material at maximum amplitude.
Dimensions (mm)	L 385 xW 305 x H 245
Net Weight	50 kg Without Sieves (Approx)

Advance Model (iControl-QS-200C)


SOP STORAGE UPTO 8 PROGRAMS

ADVANCE MODE PROGRAM OPTION



10mm
20µm

NEW



Dry Sieving	Standard
Wet Sieving	Optional Kit
Mode of Operation	3 modes (Continuous, Intermittent & Advance)
Intermittent Operation	At intervals upto 0.5 seconds
Advance Mode	Programmable time for Continuous & Intermittent (Together) in this mode
Capacity	Upto 8 Sieves of 50mm height Upto 16sieves of 25 mm height
Shake Time	Programmable from 1 mm to 99mm
Power Level	Programmable from 1 to 20
Amplitude	0.1mm to 2.5mm
Display	20 x 4 character alphanumeric LCD
Noise Level	Less than 61 db without sieves at maximum amplitude. Less than 71 db with sieves and material at maximum amplitude.
Dimensions (mm)	L 385 xW 305 x H 245
Net Weight	50 kg Without Sieves (Approx)

ARROW PHARMATECH

C-004, Vijaya CHS., D.V. Marg, Vitawa, Thane (W) - 400 605. India

Email : arrowpharmatech@gmail.com / Contact No. 8169231891 / 7400063721

URL: www.arrowpharmatech.com E-mail: arrowpharmatech@gmail.com, kumar.m@arrowpharmatech.com



ARROW PHARMATECH

ARROW PHARMATECH
EXPERTISE . SERVICE . INTEGRITY


Electromagnetic Sieve Shakers

Primary Model (QS-200B)

THREE DIMENSIONAL
SIEVING MOTION




10mm
20µm



Dry Sieving	Standard
Wet Sieving	Optional Kit
Mode of Operation	2 modes (Continuous, Intermittent)
Intermittent Operation	At intervals upto 0.5 seconds
Capacity	Upto 8 Sieves of 50mm height Upto 16sieves of 25 mm height
Shake Time	Programmable from 1 mm to 99mm
Power Level	Programmable from 1 to 20
Amplitude	0.1mm to 2.5mm
Display	16 x 2 character alphanumeric LCD
Noise Level	Less than 61 db without sieves at maximum amplitude. Less than 71 db with sieves and material at maximum amplitude.
Dimensions (mm)	L 385 xW 305 x H 245
Net Weight	50 kg Without Sieves (Approx)

Test Sieve



ISO NORMAL APERTURE SUPPLEMENTARY SIZE R40/3	(MESH SIZE) ASTM E11-01	BSS
6.0mm		
4.75mm		
4.00mm	No.5	No.4
3.35mm	No.6	No.5
2.80mm	No.7	No.6
2.36mm	No.8	No.7
2.00mm	No.10	No.8
1.70mm	No.12	No.10
1.40mm	No.14	No.12
1.18mm	No.16	No.14
1.00mm	No.18	No.16
850µm	No.20	No.18
710µm	No.25	No.22
600µm	No.30	No.25
500µm	No.35	No.30
425µm	No.40	No.36
355µm	No.45	No.44
300µm	No.50	No.52
250µm	No.60	No.60
212µm	No.70	No.72
180µm	No.80	No.85
150µm	No.100	No.100
125µm	No.120	No.120
106µm	No.140	No.150
90µm	No.170	No.170
75µm	No.200	No.200
63µm	No.230	No.240
53µm	No.270	No.300
45µm	No.325	No.350
38µm	No.400	
32µm	No.450	
25µm	No.500	
20µm	No.635	

ARROW PHARMATECH

C-004, Vijaya CHS., D.V. Marg, Vitawa, Thane (W) - 400 605. India

Email : arrowpharmatech@gmail.com / Contact No. 8169231891 / 7400063721

URL: www.arrowpharmatech.com E-mail: arrowpharmatech@gmail.com, kumar.m@arrowpharmatech.com