




ARROW PHARMATECH
EXPERTISE . SERVICE . INTEGRITY

ARROW PHARMATECH

Laboratory Glassware Washer



Aurora-F2 Open the door view

Dimension (H*W*D)	995*930*765mm
Number of cleaning layers	2 layers
Chamber volume	198L
Circulation pump flow rate	0-600L/min adjustable
Heating Power	4kw/9kw
Basket identification system	Standard
Installation method	Free standing
Washing capacity (For example)	①25ml Volumetric flask 144 seats ②100ml volumetric flask 84 seats ③Sample vials 476 seats ④Pipettes 238 seats ⑤Petri dishes 168 seats

ARROW PHARMATECH

C-004, Vijaya CHS., D.V. Marg, Vitawa, Thane (W) - 400 605. India

Email : arrowpharmatech@gmail.com / Contact No. 8169231891 / 7400063721

URL: www.arrowpharmatech.com E-mail: arrowpharmatech@gmail.com, kumar.m@arrowpharmatech.com



ARROW PHARMATECH
EXPERTISE . SERVICE . INTEGRITY

ARROW PHARMATECH

Special Cleaning Baskets


Moment / Glory / Aurora / Flash Washing Frames



FA-Z01

Upper level basket frame

With two connectors, Used to connect 2 injector inserts
Auto self-sealing docking valve
H140, W536, D562mm



FA-Z03

Upper level basket frame

TO load shelf
Height adjustable
Built-in spray arm
H183, W535, D552mm



FA-Z05

Upper level basket frame

TO load shelf
Height adjustable
Built-in spray arm
H152, W536, D567mm



FA-Z02

Lower level basket frame

With two connectors, Used to connect 2 injector inserts
Auto self-sealing docking valve
H158, W528, D566mm



FA-Z04

Lower level basket frame

To load tray and shelf
H79, W529, D546mm



FA-Z06

Lower level basket frame

Double level
Built-in spray arm
H226, W529, D569mm



FA-Z11

Full frame injection module

38 pipettes can be loaded in three layers
10 10pcs 10-100ml pipette, H550MM
14pcs 10-25ml pipettes, H500MM
14pcs 1-10ml pipettes, H500MM
H373, W528, D558mm

ARROW PHARMATECH

C-004, Vijaya CHS., D.V. Marg, Vitawa, Thane (W) - 400 605. India

Email : arrowpharmatech@gmail.com / Contact No. 8169231891 / 7400063721

URL: www.arrowpharmatech.com E-mail: arrowpharmatech@gmail.com, kumar.m@arrowpharmatech.com



ARROW PHARMATECH
EXPERTISE . SERVICE . INTEGRITY

ARROW PHARMATECH

Special Cleaning Baskets

Moment / Glory / Aurora / Flash Modules/Baskets



FA-K08

Injection Module 8 injections

For 500-2000ml Erlenmeyer flasks, volumetric flask, measuring cylinder and so on
Injector nozzle: $\Phi 6 \times H220$ mm
H253, W140, D488 mm



FA-K11

Injection Module 10 injections

For 500-2000ml Erlenmeyer flasks, volumetric flask, measuring cylinder and so on
Injector nozzle: $\Phi 6 \times H220$ mm
H253, W140, D488 mm



FA-K21

Injection Module 21 injections

For Erlenmeyer flasks, volumetric flask, measuring cylinder and so on
Injector nozzle: $\Phi 4 \times H160$ mm
H183, W192, D488 mm



FA-K31

Injection Module 36 injections

For Erlenmeyer flasks, volumetric flask, measuring cylinder and so on
Injector nozzle: $\Phi 4 \times H160$ mm
H164, W204, D511 mm



FA-K10

Injection Module 10 injections

For 500-2000ml Erlenmeyer flasks, volumetric flask, measuring cylinder and so on
Injector nozzle: $\Phi 6 \times H220$ mm
H253, W140, D488 mm



FA-K20

Injection Module 21 injections

For Erlenmeyer flasks, volumetric flask, measuring cylinder and so on
Injector nozzle: $\Phi 4 \times H160$ mm
H253, W192, D488 mm



FA-K22

Injection Module 21 injections

For Erlenmeyer flasks, volumetric flask, measuring cylinder and so on
Injector nozzle: $\Phi 2.5 \times H110$ mm
 $\Phi 4 \times H160$ mm
 $\Phi 6 \times H220$ mm
H253, W192, D488 mm



FA-K30

Injection Module 36 injections

For Erlenmeyer flasks, volumetric flask, measuring cylinder and so on
Injector nozzle: $\Phi 2.5 \times H110$ mm
H114, W204, D511 mm



FA-K80

Injection Module 84 injections


For 84pcs glass tubes, centrifuge tube, test tubes and so on.
Injector nozzle: $\Phi 2.5 \times H70$ mm
H115, W190, D503mm



FA-K90

Injection Module 119 injections

Can load 119pcs glass tubes, centrifuge tube, test tubes, sampling vials and so on.
Injector nozzle: $\Phi 2.5 \times H45$ mm
H114, W190, D490mm



FA-K90/1

Injection Module 119 injections

Can load 119pcs glass tubes, centrifuge tube, colorimetric tubes, test tubes and so on.
Injector nozzle: $\Phi 2.5 \times H90$ mm
H190, W190, D490mm



FA-K91

Injection Module 119 injections

Can load 119pcs pipettes, Pipettes' maximum length can be 460mm
H189, W190, D485mm

ARROW PHARMATECH

C-004, Vijaya CHS., D.V. Marg, Vitawa, Thane (W) - 400 605. India

Email : arrowpharmatech@gmail.com / Contact No. 8169231891 / 7400063721

URL: www.arrowpharmatech.com E-mail: arrowpharmatech@gmail.com, kumar.m@arrowpharmatech.com



ARROW PHARMATECH
EXPERTISE . SERVICE . INTEGRITY

ARROW PHARMATECH

Special Cleaning Baskets

Moment / Glory / Aurora / Flash Modules/Baskets



FA-M36

Injection Module 36 injections

Can load 28pcs pipettes, 8pcs Erlenmeyer flasks, volumetric flask, measuring cylinder and so on
Injector nozzle: $\Phi 6^*H220$ mm
H255, W190, D493mm



FA-M65

Injection Module 65 injections

Can load 28pcs pipettes, 28pcs sampling vials, 9pcs Erlenmeyer flasks, volumetric flask, measuring cylinder and so on
Injector nozzle: $\Phi 2.5^*H110$ mm; $\Phi 4^*H160$ mm; $\Phi 6^*H220$ mm
H265, W190, D493mm



FA-M112

Injection Module 112 injections

Can load 56pcs pipettes, 56pcs test tubes and so on
Injector nozzle: $\Phi 2.5^*H110$ mm;
H190, W190, D494mm



T-101

Full frame 56

Applied to clean microbiological laboratory petri dish
56 brackets, 70mm in height
H145, W504, D434mm



T-102

Full frame 56

Only used together with TXF-101
Applied to clean microbiological laboratory petri dish
56 brackets, 70mm in height
H103, W507, D414mm



T-201

Basket

28 spring clips are used to hold the vessel
H116, W220, D410mm



T-202

Basket


28 spring clips are used to hold the vessel
H186, W220, D445mm



T-203

Basket

21 spring clips are used to hold the vessel
H158, W230, D460mm



T-204/1

dividing

T-204/2

Basket

For funnel, beaker, jar, etc
H120, W224, D434mm

G-401

Cover net

for fixing delicate wash items

* for T-204



T-401/402/403/ 404

Basket

H100, W218, D218mm The maximum is 12*75mm
H127, W218, D218mm The maximum is 12*105mm
H187, W218, D218mm The maximum is 12*165mm
H230, W218, D218mm The maximum is 12*200mm



G-201/G-202

Cover net

Basket covering for sample bottle to prevent the glassware rushed out .

G-201 is suitable for: FA-K90
G-202 is suitable for: T-101/T-102
H9, W215, D455 mm



B-101

Holder

Professional holder, Makes the bottles more stable

Z-101

Pipette cleaning sleeve

For pipettes
 $\Phi 8^*H116$ mm

ARROW PHARMATECH

C-004, Vijaya CHS., D.V. Marg, Vitawa, Thane (W) - 400 605. India

Email : arrowpharmatech@gmail.com / Contact No. 8169231891 / 7400063721

URL: www.arrowpharmatech.com E-mail: arrowpharmatech@gmail.com, kumar.m@arrowpharmatech.com



ARROW PHARMATECH
EXPERTISE . SERVICE . INTEGRITY

ARROW PHARMATECH

Special Cleaning Baskets

Moment / Glory / Aurora / Flash Modules/Baskets



R-201

Combination structure basket frame

Double level, fit to injection and non injection basket.
Heating air and water pumped into basket through connection system.



R-203

Combination structure basket frame

Three levels, including two spray arms
Connecting drying system With Basket recognition system.



R-204

Combination structure basket frame

Four levels, including three spray arms
Connecting drying system
With Basket recognition system.



R-205

Combination structure basket frame


Five levels, including four spray arms
Connecting drying system
With Basket recognition system.



SX-901

Injection basket module

Suitable for volumetric flask, colorimeter tubes, Erlenmeyer flasks, colorimeter flask, flat bottom flask etc.



SX-902

Injection basket module

Suitable for centrifuge tube, sampling vials, test tube etc.



DZ-901

Injection basket module

Suitable for volumetric flask, Erlenmeyer flasks, colorimeter flask, flat bottom flask etc.



DZ-902

Injection basket module

Suitable for Pipettes, the maximum height is 580mm, can wash 116 injections



T-480

Trolley

Connected to machine, for loading and supporting baskets

ARROW PHARMATECH

C-004, Vijaya CHS., D.V. Marg, Vitawa, Thane (W) - 400 605. India

Email : arrowpharmatech@gmail.com / Contact No. 8169231891 / 7400063721

URL: www.arrowpharmatech.com E-mail: arrowpharmatech@gmail.com, kumar.m@arrowpharmatech.com