



**ARROW PHARMATECH**  
EXPERTISE . SERVICE . INTEGRITY

# ARROW PHARMATECH

## PH ELECTRODES CATALOGUE

*Why use ordinary when you have the Best in Class.....*



### AGL-42 Glass Body 5 PT

<b>MODEL</b>	<b>AGL -42</b>	<b>HIGHLIGHTS</b> Suitable for 5 point calibration  Get Stable Reading to 0.02 ph in less than 30 Sec.  Results are 3 to 5 times more Precise.  Available with Teflon protective cap
<b>DESCRIPTION</b>	<b>Glass Electrode</b>	
<b>RANGE</b>	<b>0 – 14 PH</b>	
<b>TEMPERATURE</b>	<b>0 -90 °C</b>	
<b>TYPE REF.</b>	<b>Ag / AgCl</b>	
<b>JUNCTION</b>	<b>Ceramic</b>	
<b>REFFILABLE</b>	<b>Yes</b>	
<b>BODY</b>	<b>Glass</b>	
<b>CONNECTOR</b>	<b>BNC</b>	
<b>SIZE (MM)</b>	<b>160 x 12</b>	
<b>CABLE</b>	<b>1 meter</b>	

### APY-41 EPOXY Electrode

<b>MODEL</b>	<b>APY- 41</b>	<b>HIGHLIGHTS</b> Suitable for 3 point calibration  Get quicker stable readings  Useful for aqua's solution samples  With built in Protective cap
<b>DESCRIPTION</b>	<b>Plastic electrode</b>	
<b>RANGE</b>	<b>0 – 14 PH</b>	
<b>TEMPERATURE</b>	<b>0 -80 °C</b>	
<b>TYPE REF.</b>	<b>Ag / AgCl</b>	
<b>JUNCTION</b>	<b>Fiber</b>	
<b>REFFILABLE</b>	<b>No</b>	
<b>BODY</b>	<b>Poly Carbonate</b>	
<b>CONNECTOR</b>	<b>BNC</b>	
<b>SIZE (MM)</b>	<b>160 x 12</b>	
<b>CABLE</b>	<b>1 meter</b>	



## ARROW PHARMATECH

C-004, Vijaya CHS., D.V. Marg, Vitawa, Thane (W) - 400 605. India

Email : [arrowpharmatech@gmail.com](mailto:arrowpharmatech@gmail.com) / Contact No. 8169231891 / 7400063721

URL: [www.arrowpharmatech.com](http://www.arrowpharmatech.com) E-mail: [arrowpharmatech@gmail.com](mailto:arrowpharmatech@gmail.com), [kumar.m@arrowpharmatech.com](mailto:kumar.m@arrowpharmatech.com)



**ARROW PHARMATECH**  
EXPERTISE . SERVICE . INTEGRITY

# ARROW PHARMATECH

## PH ELECTRODES CATALOGUE

*Why use ordinary when you have the Best in Class.....*

### APC-22 Glass Body 3 PT



<b>MODEL</b>	APC-22
<b>DESCRIPTION</b>	Glass Electrode
<b>RANGE</b>	0 – 14 PH
<b>TEMPERATURE</b>	0 -90 °C
<b>TYPE REF.</b>	Ag / AgCl
<b>JUNCTION</b>	Ceramic
<b>REFFILABLE</b>	Yes
<b>BODY</b>	Glass
<b>CONNECTOR</b>	BNC
<b>SIZE (MM)</b>	160 x12
<b>CABLE</b>	1 meter

#### HIGHLIGHTS

Suitable for 3 point calibration

Get Stable Reading to 0.05 ph in less than 45 sec.

Cost effective and efficient results.

Available with Teflon protective cap

### PH-2100 Ultem Body

<b>Diameter</b>	12mm
<b>Length</b>	150mm
<b>Reference</b>	Sealed
<b>Junction Number</b>	Double (2)
<b>Junction Material</b>	HDPE & Pellon
<b>Reference Solution</b>	3.5 M KCl/AgCl /KNO <sub>3</sub> (gel)
<b>Glass Shape</b>	Flat
<b>Body Material</b>	Ultem
<b>Temperature Range</b>	0-80 degrees C
<b>Strain Relief</b>	No
<b>BNC Boot</b>	Yes



### PH-2500 Flat Surface

<b>Diameter</b>	12mm
<b>Length</b>	150mm
<b>Reference</b>	Sealed
<b>Junction Number</b>	Double (2)
<b>Junction Material</b>	HDPE & Pellon
<b>Reference Solution</b>	3.5 M KCl/AgCl /KNO <sub>3</sub> (gel)
<b>Glass Shape</b>	Flat
<b>Body Material</b>	Glass
<b>Temperature Range</b>	0-80 degrees C
<b>Strain Relief</b>	No
<b>BNC Boot</b>	Yes



## ARROW PHARMATECH

C-004, Vijaya CHS., D.V. Marg, Vitawa, Thane (W) - 400 605. India

Email : arrowpharmatech@gmail.com / Contact No. 8169231891 / 7400063721

URL: www.arrowpharmatech.com E-mail: arrowpharmatech@gmail.com, kumar.m@arrowpharmatech.com